

Temporary Tattoo



1. Print Tattoo Sheet in Mirror Image mode.
 2. Remove Release Liner from Glue Film.
 3. Place printed Tattoo Sheet Image down on the Glue Film. Rub firmly with Cloth or run through cold Laminator.
 4. Trim Tattoo Transfer to final shape and remove clear Film from Glue.
 5. Place in desired Position on Skin and cover with wet Cloth i.e. Towel.
 6. After approx. 20 Seconds slide the BasePaper sideways and remove. Softly wipe off the excess release residue.
- Done.